



Measuring Capacity


Name _____

Score _____


MC:34

Use the readings given below to shade each cylinder.


1) 84 mL




2) 420 mL




3) 99 mL




4) 3.6 mL




5) 325 mL




6) 7.9 mL



7) 160 mL



8) 12 mL





Measuring Capacity

Name _____


Score _____

Answer key


MC:34

Use the readings given below to shade each cylinder.


1) 84 mL




2) 420 mL




3) 99 mL




4) 3.6 mL




5) 325 mL



6) 7.9 mL



7) 160 mL



8) 12 mL

