



Measuring Capacity


Name _____

Score _____


MC:35

Use the readings given below to shade each cylinder.


1) 290 mL




2) 4.3 mL




3) 500 mL




4) 76 mL




5) 3.1 mL




6) 52 mL



7) 48 mL



8) 365 mL





Measuring Capacity

Name _____


Score _____

Answer key


MC:35

Use the readings given below to shade each cylinder.


1) 290 mL




2) 4.3 mL




3) 500 mL




4) 76 mL




5) 3.1 mL



6) 52 mL



7) 48 mL



8) 365 mL

